



# Microsoft Rapid Screening Solution

V4 (2.0.0.7) Deployment Guide

April 26, 2021



## VERSION 4 (2.0.0.7)

This release we are introducing home screening. Employees will be able to take advantage of a step-by-step approach guiding them from the intake of a sample through the input of the result. In the solution settings we are making it easy to enable home screening and input the desired consent verbiage.

Read on to find out the feature details!

# Version 4 (2.0.0.7) - Features

## Employee App

- Home screening new menu
- Employee home screening consent
- Employee home screening attestation
- Home screening steps by steps videos
- Home screening timer
- Home screening result view
- Home screening result photo capability

## Solution Settings

- Allow home screening
- Need to attest for home screening
- Photo mandatory for home screening
- Home screening consent


## Environment Variables

- Home screen wait time in minutes
- Home screen time out in minutes


## Dataverse

- Home screening table to store steps, images and video links

# Home Screening - User Journey



# Home Screening - Demo




## Notes:

1. The minimum timer was set to a few seconds for demo purposes.
2. The picture was taken with a pen for demo purposes, and it should be a test in reality.

# Employee App




Home screening new menu



# Employee App

Employee must consent to home screening.

Through the Employee profile



Welcome Karine!

Full Name

Karine Nadeau

Email

admin@M365x216739.onmicrosoft.com

Personal Email

Phone Number

514-249-7060

\* Age Range

26-35

\* Privacy Consent

☒

\* Consent to home screening

☒

\* Workplace Location

Home screening - Quebec

Confirm

or

By clicking on HOME SCREENING

Do you agree with the following statements?

Home screening consent


I AGREE

I DISAGREE


# Employee App

## Home screening steps-by-steps videos

Note: Videos can be made by your organization and be stored in either Microsoft Stream or Azure Media Services. If you prefer to use videos from an external source you can select YouTube. More details can be found in the "V4 Deployment Guide".



## Example using YouTube








# Employee App

Employee must attest to home screening

Note: Replace "Attestation verbiage by your organization attestation wording".



Welcome Karine!

Select Test Type

Abbott Panbio




1. Prep



2. Swab



3. Swirl






4. Drop



5. Result

1

Start Screening



Welcome Karine!

2

Do you agree with the following statements?

Attestation verbiage

All future appointments will cancel if you disagree

I AGREE

I DISAGREE


# Employee App

1. Home screening timer  
(result input minimum and maximum can be set).

2. Result view


*ie. It takes minimum 15 minutes for the test to provide a result which is valid for a maximum of 2 hours.*

Note: View "V4 Deployment Guide" to find out how to set it up.



The screenshot shows the 'Home Screening Timer' interface. At the top, there's a header with a home icon, the 'CREATIVE3 DESTRUCTION Rapid Screening Consortium' logo, and a user profile icon. Below the header, it says 'Welcome Karine!'. The main content area is divided into three sections: 1. 'Take your sample', 2. 'Start timer', and 3. 'Input your result'. In the 'Start timer' section, there's a stopwatch icon showing '15:00' and a 'Start timer' button. A red circle with the number '1' is next to the stopwatch. In the 'Input your result' section, there's a 'Result' dropdown menu with a flask icon and a 'Photo of result' section with a camera icon and a 'Take photo' button. At the bottom, there's a 'REVIEW & SAVE' button.

Click on VIEW COVID SCREENING RESULTS to see this view






The screenshot shows the 'Your Bookings' interface. At the top, there's a header with a home icon, the 'CREATIVE3 DESTRUCTION Rapid Screening Consortium' logo, and a user profile icon. Below the header, it says 'Welcome Karine!' and 'Your Bookings'. The main content area shows a list of bookings. The first booking is 'Home Screening' on '4/23/2021 2:46 PM' with an 'Input result' button. The second booking is 'Timed out' 'Home Screening' on '4/22/2021 12:25 PM'. A red circle with the number '2' is next to the 'Timed out' booking. An orange line connects the 'Input result' button to the 'Timed out' booking, indicating that the user cannot input a result for a timed-out booking.

Within the threshold so the user can still input their result

Outside of the threshold so the user can't input their result and must start another screening

# Employee App


Home screening result  
photo capability



Welcome Karine!  
Home Screening

1. Take your sample

2. Start timer



00:00

Wait

Time is up

3. Input your result

Result





 

Photo of result

 **1** Take photo

REVIEW & SAVE

**2**



SAVE & GO BACK

RETAKE

# Solution Settings

## 1. Allow home screening

Do you want to allow employees to do home screening yes/no?

## 2. Need to attest for home screening

Must the employee performing home screening attest their symptoms?

## 3. Photo mandatory for home screening

Must the employee performing home screening provide a photo of the test itself?

## 4. Home screening consent

What is the home screening consent verbiage?

Registrations	No
Allow Share Guest Registration	No
1 Allow Self Screening	Yes
Mask Employee Info In Pass	No
Mask Employee Info In Screening App	No
Hide Phone Number For Employees	No
Hide Personal Email For Employees	No
2 Need To Attest For SelfScreening	Yes
3 Photo Mandatory For Self Screening	Yes

4 Self Screening Consent

Self-Screening Consent.....

Opt Out Of Program Message

\* You have not provided your consent to the collection, use and sharing of your personal data. In order to continue using the Rapid Screening Application please change your Privacy Consent answer to Yes. Otherwise please exit the application.

Email Body For Result Email

\* Thank you for participating to the Employee Screening Pilot Program. Please be advised that your screening result is now available. Please review the attached document for your result.

# Dataverse

Home screening table to store steps, images and video links

Note: View "V4 Deployment Guide" to find out how to set it up.



Display name	Name	Created	Version	Managed external dependencies	Solution check
Consortium Member Global	ConsortiumMember...	4/22/2021	2.0.0.7		Hasn't been run
Consortium Member	ConsortiumMember	2/23/2021	1.1.20200205.3		Hasn't been run
Return to the Workplace	returntoworkplace	2/18/2021	1.4.282.1		Checked by Micro...
Return to the Workplace Solutions	returntoworkplaceso...	2/18/2021	1.4.282.1		Checked by Micro...
Power Apps Checker Base	msdyn_PowerAppsC...	2/13/2021	1.2.0.34		Checked by publ...
Power Apps Checker	msdyn_PowerAppsC...	2/13/2021	1.2.0.34		Checked by publ...
Contextual Help Base	msdyn_ContextualH...	2/13/2021	1.0.0.12		Checked by publ...
Contextual Help	msdyn_ContextualH...	2/13/2021	1.0.0.12		Checked by publ...
Common Data Services Default Solution	Cr588b4	2/13/2021	1.0.0.0		Hasn't been run
Default Solution	Default	2/13/2021	1.0		Not supported f...

# Environment Variables

1. Home screen wait time in minutes
2. Home screen time out in minutes

Note: View "V4 Deployment Guide" to find out how to set it up.

1

### Edit SelfScreen Wait Time In... ✕

Use this variable to store information about an app or flow. Its values can be updated as it moves to different environments. [Learn more](#)

**Display name \***

SelfScreen Wait Time In Minutes

**Name \* ⓘ**

cdl\_SelfScreenMinWaitTime

**Description**

Minimum Self Screen Wait time in minutes

**Data Type \***

Text

**Default Value ⓘ**

15

**Current Value**

Override the default value by setting the current value for your environment.

.10 ...

ⓘ Remove this value before exporting if it shouldn't be used in other environments

Save Cancel

2

### Edit SelfScreen Time Out In ... ✕

Use this variable to store information about an app or flow. Its values can be updated as it moves to different environments. [Learn more](#)

**Display name \***

SelfScreen Time Out In Minutes

**Name \* ⓘ**

cdl\_SelfScreenTimeOutInMinutes

**Description**

The number of minutes the end user needs t...

**Data Type \***

Text

**Default Value ⓘ**

120

**Current Value**

Override the default value by setting the current value for your environment.

2 ...

ⓘ Remove this value before exporting if it shouldn't be used in other environments

Save Cancel

# Feedback

We are always looking to improve our documentation. If you have feedback and/or suggestions, please email us at:  
[RapidScreeningSolution@service.microsoft.com](mailto:RapidScreeningSolution@service.microsoft.com)

